

## ATM-INNOVATION:

### DOSING SYSTEM FOR HOT MOUNTING MATERIALS

High demands of samples and therefore consumption of a lot of mounting material recommend the ATM-dosing system.

The swivel arm is either fixable at the media rail of system lab cabinets or on the wall for use with table devices.

The dosing unit is equipped with an adapter for 5-Liter-buckets and can be taken out completely.

The dosing flask portions the mounting material adjustable from 10–30 cm<sup>3</sup>. The dosing head can be directly attached on the feeder chute of the mounting press. This minimizes soiling of the mounting mould by mounting material dust.

A double dosing system is available for using different mounting materials.



### SWIVEL SYSTEM (select equipment)

» Fixing on media rail of cabinet

Order No.: A0707000

### SWIVEL SYSTEM FOR WALL FIXING

» For cabinets without media rail

Order No.: A0707001

### DOSING SYSTEM FOR HOT MOUNTING MATERIALS

» For portioning hot mounting material (adjustable 10–30 cm<sup>3</sup>)

» Dosing unit can be taken out completely

» Adapter system for 5-Liter-bucket

Order No.: A5800439

### ADAPTER-SYSTEM FOR 5L-BUCKET

